




carmo
Your partner in plastic components

09-828 PVC/PE T-Tap Cross Valve

Simple one hand operation. Drains ¼ US-gallon in 21 seconds / 1 litre in 22 seconds



Nominal measurements in mm

Simple one hand operation. Drains ¼ US-gallon in 21 seconds / 1 litre in 22 seconds

T-Tap cross valves are typically used as bottom valves in medical and Urinebags.

Materials: PVC+PEHD

Std. packaging: 1 x 2.000

Interested in other colors or materials? Please call for quote...



Carmo A/S, Hoejvangen 19, DK-3060 Espergaerde, Denmark
T: +45 4912 2100, F: +45 4912 2199, E: denmark@carmo.dk
www.carmo.dk

More info <http://www.carmo.dk/en/medical-components/drain-valves/t-tap-valve/09-828/>

All geometrical dimensions and functional information are indicative. (Unless tolerances are specified). Carmo is not responsible for misprints and technical changes.
2024-06-17 18:00:29 UTC