




**carmo**  
Your partner in plastic components

## 09-828 PVC/PE T-Tap Cross Valve

Simple one hand operation. Drains ¼ US-gallon in 21 seconds / 1 litre in 22 seconds



Nominal measurements in mm

Simple one hand operation. Drains ¼ US-gallon in 21 seconds / 1 litre in 22 seconds

T-Tap cross valves are typically used as bottom valves in medical and Urinebags.

Materials: PVC+PEHD

Std. packaging: 1 x 2.000

Interested in other colors or materials? Please call for quote...



Carmo A/S, Hoejvangen 19, DK-3060 Espergaerde, Denmark  
T: +45 4912 2100, F: +45 4912 2199, E: [denmark@carmo.dk](mailto:denmark@carmo.dk)  
[www.carmo.dk](http://www.carmo.dk)

More info <http://www.carmo.dk/en/medical-components/drain-valves/t-tap-valve/09-828/>

All geometrical dimensions and functional information are indicative. (Unless tolerances are specified). Carmo is not responsible for misprints and technical changes.  
2022-11-27 06:33:28 UTC